


Canine Atrial RNC00-SM

**Simucore Model Based Upon: Ramirez, Nattel,
Courtemanche Model of Canine Atrial Cardiac Action
Potentials, 2000; v. 1.4**

Table of contents

1 Description.....	2
2 References.....	2
3 Ordering.....	3

1. Description



Example action potential (BCL = 400 ms).

This model simulates canine atrial action potentials. The model's formulation is based on Luo-Rudy II guinea-pig ventricular model with currents modified to mimic the configuration of canine atrial action potential.

Abstract excerpt: *"To obtain insights into the ionic mechanisms governing canine atrial AP properties, we incorporated formulations of $K(+)$, $Na(+)$, $Ca(2+)$, and $Cl(-)$ currents, based on measurements in canine atrial myocytes, into a mathematical model of the AP."*

2. References

- Ramirez RJ, Nattel S, Courtemanche M.
Mathematical analysis of canine atrial action potentials: rate, regional factors, and electrical remodeling.

Am J Physiol Heart Circ Physiol. 2000 Oct;279(4):H1767-85.
PMID: [11009464](#)

3. Ordering

- [Order this model](#) or [request further information](#).